

# Understanding the Genes in SC/SI Individuals of *Physalis acutifolia*

## Transcriptome of *Physalis acutifolia*

Purpose: To have a published transcriptome (which currently doesn't exist), to be able to untangle pleiotropic genes (Mulato-Brito, 2007), and sort out housekeeping genes.

### Approach:

1. RNA extraction of SC individual of leaf, fruit, flower, root, and bud.
2. Sequence this in the Nanopore with a coverage of 80x.
  - a. This would require 1 flow cell for all 6 samples.

Genome size ~3.5 Gb.

Estimate transcriptome to be 1% of that ( $0.035 \text{ GB} \times 80 = 2.8 \text{ Gb}$ ) for one tissue. 6 different tissues at ~16.8 Gb. Nanopore ranges from 10-20 Gb.

1. De Novo Assembly
  - a. Steps below
2. Annotation of transcriptome
3. Potentially publish at: ???

## Differentially Expressed Genes (DEG) between SC/SI *Physalis acutifolia*

Purpose: To understand the genes that regulate SI within this species.

### Approach:

1. RNA extraction of unpollinated styles.
  - a. Practice with 2 samples. 1 SC and 1 SI of one individual with 5 styles pooled.
    - i. Optimize to obtain good yield.
  - b. Do 2 samples with pollinated SC/SI and compare with unpollinated.
    - i. After talking with Pat she said there wasn't much difference with pollinated and unpollinated styles.

**\*\*Reevaluate\*\*** methods

2. Run Scotty - Power Analysis for RNA Seq Experiments with the samples to reevaluated number of samples needed.

- a. Continue with 3 populations (1 SC and 2 SI). 3 individual plants with 5 styles pooled together.

**\*Evaluate** data and revise methods if needed.

- a. Finish with the next 7 populations -- 3 individual plants with 5 styles pooled together.

Total Extractions: 32 extractions

3. De Novo Assembly
4. Determine DEG
5. qPCR on DEG
5. Potentially publish results at: AJB, Genetics

## Sequencing Protocol

1. Extract RNA: Qiagen Plant RNEasy Kit
  - a. Qiagen DNA Cleanup during extraction
2. Check Quality:
  - a. Nanodrop for quality
  - b. Bioanalyzer
3. Clean-up RNA: NEBNext Ploy(A) mRNA Magnetic Isolation Module
4. cDNA Synthesis: Native Barcoding Expansion 1-12 (PCR-free) and add in oligo(dT)/Random hexamer primers
5. Library prep and sequencing: Oxford Nanopore protocols and sequence using MinION SpotON Flow Cell

## Computational Protocol

1. Sequence Processing: Basecalling using Albacore (Nanopore tech.)
2. Demultiplexing using porechop
3. DeNovo Assembling: BWA, SAMtools, SOAPdenovo-Trans
4. Determine read depth: External RNA Controls Consortium RNA Spike-In Mix
5. Transcript analysis: PANTHER Overrepresentation with Gene Ontology database
6. DEseq R to determine DEG

- 
1. With DEG do qPCR

## Cost Breakdown

NEBNext Poly(A) mRNA magnetic isolation (24 samples)  
 $\$67 \times 2 = \$134$

Nanopore Kit  
 $\$1000.00 \times 2 = \$2000.00$

SuperScript IV (50 reactions)  
 $\$488.00$

PerfeCTa SYBR® Green SuperMix  
 $\$194.60$

Gene Specific Primers (~10 genes including control primers)  
 $\$150.00$

Shipping Cost  
 $\$150$

*Total: \$3116.60*